

# CENTRIC CF 220

## CLINICAL CENTRIFUGES



### MAIN FEATURES:

- > Ergonomic design
- > User friendly
- > Stainless steel chamber
- > Quick parameter settings
- > Large clearly visible parameters

Ergonomic and user friendly clinical centrifuge **Centric CF 220** is instrument for clinical use and comes with eight place angle rotor for tubes from 3 ml to 15 ml.

Large illuminated LCD display offers great visibility of set parameters.

Centrifuge **Centric CF 220** is simple model, easy to operate. Its stainless steel chamber offers easy cleaning of the centrifuge's interior.

| TECHNICAL DATA                   | CENTRIC CF 220                                                    |
|----------------------------------|-------------------------------------------------------------------|
| Max. rotation speed              | 200 to 6.000 rpm in steps of 100 rpm                              |
| Max. centrifugal force           | 100 do 4.180 x g in steps of 100 x g                              |
| Max. capacity                    | 8 x 15 ml                                                         |
| Time                             | 10s to 99:50min, in steps of 10s, continuous run at HOLD function |
| Programs                         | 100                                                               |
| Acceleration & Deceleration      | 10 / 10 levels                                                    |
| Noise level (at max. speed)      | ≤ 59 dBA                                                          |
| Fast spin with preselected speed | Yes                                                               |
| Permitted ambient conditions     | Temperature: 2°C to 35°C,<br>Humidity: 75%, non-condensing        |
| Operation control                | Microprocessor                                                    |
| Motor type                       | Brushless induction motor, maintenance free                       |
| Dimensions (W x D x H)           | 280 x 370 x 226 mm                                                |
| Weight                           | 13,9 kg                                                           |

| ORDERING INFORMATION           | Order ref. No. | DOMEL code   |
|--------------------------------|----------------|--------------|
| <b>CENTRIC CF 220</b>          |                |              |
| 230V / 50 - 60 Hz              | 442742         | 232.005      |
| <b>CENTRIC CF 220</b>          |                |              |
| 230V / 50 - 60 Hz              | 466524         | 232.005/UK   |
| <b>CENTRIC CF 220 120V</b>     |                |              |
| 120V / 50 - 60 Hz              | 445418         | 1232.005     |
| <b>CENTRIC CF 220 IVD</b>      |                |              |
| 230V / 50 - 60 Hz              | 465561         | 232.005/IVD  |
| <b>CENTRIC CF 220 IVD 120V</b> |                |              |
| 120V / 50 - 60 Hz              | 465562         | 1232.005/IVD |

We reserve the right to alter specification details etc. without prior notice or liability!

V:11-21

# DOMEL

DOMEL, d.o.o., Otoki 21, 4228 Železniki  
Domel LABS, Na Plavžu 79, 4228 Železniki

Phone: +386 4 5117 525  
E-mail: [info.labs@domel.com](mailto:info.labs@domel.com)  
WEB: <http://labs.domel.com>



Domel LABS